



Safe and reliable optical processing

Baskets for optics

Personalized protective baskets for the organized and safe packaging of scopes


+ Durable
+ Safe

+ Economical
+ Versatile




The FAMI processing basket replaces perforated trays and wire baskets

Secure closure



Adjustable silicone supports



Packaging System



FAMI Sterile Goods Basket



1 Associated with the **Flexible Packaging** it can be accommodated in the **sterile goods basket**





2 Associated with the container forms a perfect **Rigid Packaging System**

DIMENSIONS





Dimensions	Codes	Specifications
290 x 80 x 50 mm	124.111	Basket for straight rigid optics up to 28cm, with silicone supports
460 x 80 x 50 mm	124.112	Basket for straight rigid optics up to 45cm, with silicone supports
670 x 80 x 50 mm	124.113	Basket for straight rigid optics up to 66cm, with silicone supports

Dimensions	Codes	Specifications
290 x 160 x 50 mm	124.114	Basket for straight rigid optics up to 28cm, camisa and trocater, with silicone supports
460 x 160 x 70 mm	124.115	Basket for straight rigid optics up to 45cm, camisa and trocater, with silicone supports
460 x 160 x 70 mm	124.116	Basket for straight rigid optics up to 45cm, and instruments, with silicone supports




Dimensions	Codes	Specifications
290 x 80 x 50 mm	124.121	Basket for 2 straight rigid optics up to 28cm, with silicone supports
460 x 80 x 50 mm	124.122	Basket for 2 straight rigid optics up to 45cm, with silicone supports

Dimensions	Codes	Specifications
540 x 250 x 70 mm	124.131	Basket for flexible optics, with silicone supports



Dimensions	Codes	Specifications
540 x 250 x 70 mm	126.312	Basket for angled rigid optics, with silicone supports
640 x 150 x 80 mm	124.141	Basket for angled rigid optics, with silicone supports

Dimensions	Codes	Specifications
290 x 80 x 50 mm	124.211	Basket for straight rigid optics up to 28cm, without silicone supports
460 x 80 x 50 mm	124.212	Basket for straight rigid optics up to 45cm, without silicone supports



AUTHORIZED DISTRIBUTOR

