

✓ **The first step towards a safer, modern and sustainable CSSD**

Processing trays/baskets


Optimizes manual and automated cleaning process, drying, sterilization and transport.

- ✚ **Resistent**
- ✚ **Safe**
- ✚ **Lightweight**
- ✚ **Cost-effective**
- ✚ **Versatile**
- ✚ **Ergonomic**

- The FAMI processing basket replaces** perforated trays and wire baskets
- Small holes** on the sides prevent instruments from coming out.
- Large holes** in the bottom and lid facilitate the circulation of disinfecting and sterilizing agents.




Rounded corners

Prevent damage to flexible packaging



Multifunction Retractable Handles

Improves ergonomics and facilitates aseptic presentation



Washer Disinfector

Packaging System



Sterile Goods Basket



Container System



1 Associated with a **flexible sterile barrier system** can be accommodated in a **sterile goods basket**

2 Associated with a rigid sterile barrier system (container) becomes the **best Packaging System**

NEW Processing
Trays/Baskets



Codes and Dimensions of Baskets - DIN International Standard

Dimension	PN	Dimension	PN	Dimension	PN
540 x 255 x 30 mm	662.130	405 x 255 x 30 mm	662.330	245 x 255 x 30 mm	662.530
540 x 255 x 50 mm	662.150	405 x 255 x 50 mm	662.350	245 x 255 x 50 mm	662.550
540 x 255 x 70 mm	662.170	405 x 255 x 70 mm	662.370	245 x 255 x 70 mm	662.570
540 x 255 x 100 mm	662.190	405 x 255 x 100 mm	662.390	245 x 255 x 100 mm	662.590

Dimension	PN	Dimension	PN	Dimension	PN
485 x 255 x 30 mm	662.230	350 x 255 x 30 mm	662.430	265 x 120 x 45 mm	662.745*
485 x 255 x 50 mm	662.250	350 x 255 x 50 mm	662.450	265 x 170 x 45 mm	662.845*
485 x 255 x 70 mm	662.270	350 x 255 x 70 mm	662.470		
485 x 255 x 100 mm	662.290	350 x 255 x 100 mm	662.490		


Codes and Dimensions of Lids**

Dimension	PN
540 x 255 mm	662.110
485 x 255 mm	662.210
405 x 255 mm	662.310
350 x 255 mm	662.410
245 x 255 mm	662.510




* non available with lid
**lids sold separately

NEW Small Parts
Basket




- ✓ Developed to protect small parts, parts and accessories, during and after processing;
- ✓ Free of wire mesh;
- ✓ Comes with lid.

Dimension	PN
80 x 50 x 40 mm	134.117
80 x 80 x 40 mm	134.120
150 x 140 x 40 mm	134.133



NEW Sterile Goods
Basket

Organizes and improves workflow




- ✓ Optimizes processing areas;
- ✓ Improves the organization and visibility of stock;
- ✓ Ergonomic and lightweight; easy to handle;
- ✓ Increases instrument validity;
- ✓ Reduces the risk of infections.

The **side opening** guarantee the removal of the sterile packaging system with **no risks for packages and professionals.**

Dimension	PN
600 x 292 x 150 mm	665.905
600 x 292 x 260 mm	665.906

AUTHORIZED DISTRIBUTOR



**Stackable
and Pluggable**

The standardization of sizes makes it possible **to stack the baskets during assembly**, in addition **to organizing loads in the autoclaves and storing them in the arsenal.**

**Ergonomically
designed**

