



Steam Cleaner

entahal
medical

EXCLUSIVE REPRESENTATIVE
FAMI

Steamer

A new **fast and efficient** cleaning technology that uses a steam jet to fluidize and drag organic and inorganic soil from **simple, complex and cannulated** medical devices.

Remove easily organic and inorganic soil strongly adhered: substances or oily waste and stickers.



Easy supply



Professional protection



Productivity and efficiency



Sustainable technology

Multi-gun with Luer-lock connection and 4 functions:

- ✓ Dry steam;
- ✓ Hot water under pressure;
- ✓ Wet steam;
- ✓ Compressed air.



Multisteam
Manual water tank filling

Autosteam
Automatic water tank filling

Steambox
Suction/Filtering station

Descriptions



Description	PN	Dimensions	Capacity
Multisteam	8800681	41 x 37 x 35 cm	5 L
Autosteam	8800701	41 x 37 x 35 cm	5 L
Steambox	8800682-2	60 x 50 x 45 cm	-

Cleaning brushes




Description	PN	Dimensions	Quantity
S-Flex Hollow Steam Brush complete	LT635-25-S	Ø 6-8 mm, 25 cm	1 unit
Refill Pack S-Flex Steam Brush	LT635-S	Ø 6-8 mm	3 units
S-Flex Hollow Steam Brush complete	LT655-25-S	Ø 10-13 mm, 25 cm	1 unit
Refill Pack S-Flex Steam Brush	LT655-00-S	Ø 10-13 mm	3 units

Cleaning brushes External surfaces



Description	PN	Dimensions	Quantity
Steam Brush	LT581-S	Ø 20 mm	1 unit

Steambox Filter



Description	PN	Quantity
Filter	8800	5 units

Face shield



Description	PN	Dimensions	Quantity
Face shield	05TD3119	-	100 pieces
Nitrile gloves	PS-3224	Size Small	1 pair
Nitrile gloves	PS-3324	Size Medium	1 pair
Nitrile gloves	PS-3424	Size Large	1 pair

Hand, forearms and arms protector

